

ADVANTAGES OF PTFE

- 100% Pure PTFE (*without any dyes or other additives*)
- Excellent biocompatibility
- Ideally suited for application of implants
- Very soft to avoid any damage
- High knot pull strength
- Free of APFO, absolutely no health hazard
- Surface structure prevents merging of skin and tissues
- Hydrophobic and oleophobic (*no blood penetration*)
- Extremely low friction coefficient (*for skin and tissue protection*)
- Little or no package memory

“PROFIMED is a non absorbable, monofilament suture made from polytetrafluoroethylene (PTFE) polymer.”

PROFIMED
 • Undyed monofilament non absorbable, PTFE suture
 • PTFE chir. Nahtmaterial, monofil., nicht resorbierbar

6/0 USP
 (0.7 metric)
 (Eur. Pharm)
1 x 45 cm

Reverse Cutting Needle
 SCHNEIDENDE NADEL (DS)

REF PR620

16 mm CIRCLE **KREIS 3/8**

medipac
 MEDICAL SUPPLIES

Industrial Area
 61 100 Kilkis - Greece
 Industriegebiet
 61 100 Kilkis Griechenland
 Manufactured In Greece
 Hergestellt in Griechenland

STERILE EO

CE 0653

PROFIMED PTFE SUTURES

PTFE Sutures - PROFIMED is an ideal suture for dental bone crafting and implant procedures where a soft monofilament suture is desirable.

PROFIMED Suture are most often used in conjunction with bone crafting, the placement of dental implants, periodontal surgery and soft tissue grafts.

PROFIMED is biologically inert and chemically non-reactive. In addition, the monofilament construction prevents bacterial wicking which is associated with braided sutures.

Unlike other synthetic monofilament sutures, the material is supple and well tolerated in the mouth.



Reverse Cutting
 Schneidende Nadel (DS)

3/8 CIRCLE
 KREIS 3/8

▼ 12mm

| METRIC | USP | SUTURE LENGTH FADEN-LÄNGE | CODE ARTIKEL |
|--------|-----|------------------------------|-----------------|
| 0.7 | 6/0 | 45cm | PR610 |

▼ 16mm

| METRIC | USP | SUTURE LENGTH FADEN-LÄNGE | CODE ARTIKEL |
|--------|-----|------------------------------|-----------------|
| 0.7 | 6/0 | 45cm | PR620 |
| 1.0 | 5/0 | 45cm | PR621 |
| 1.5 | 4/0 | 45cm | PR622 |
| 2.0 | 3/0 | 45cm | PR623 |

▼ 19mm

| METRIC | USP | SUTURE LENGTH FADEN-LÄNGE | CODE ARTIKEL |
|--------|-----|------------------------------|-----------------|
| 0.7 | 6/0 | 45cm | PR630 |
| 1.0 | 5/0 | 45cm | PR631 |
| 1.5 | 4/0 | 45cm | PR632 |
| 2.0 | 3/0 | 45cm | PR633 |



Round Bodied - Diamond Point
 Schneidende rundkörper
 Nadel (DRT)

3/8 CIRCLE
 KREIS 3/8

⊙ 18mm

| METRIC | USP | SUTURE LENGTH FADEN-LÄNGE | CODE ARTIKEL |
|--------|-----|------------------------------|-----------------|
| 1.0 | 5/0 | 45cm | PR231T |
| 1.5 | 4/0 | 45cm | PR232T |



Reverse Cutting
 Schneidende Nadel (HS)

1/2 CIRCLE
 KREIS 1/2

▼ 18mm

| METRIC | USP | SUTURE LENGTH FADEN-LÄNGE | CODE ARTIKEL |
|--------|-----|------------------------------|-----------------|
| 1.5 | 4/0 | 45cm | PR522 |

▼ 22mm

| METRIC | USP | SUTURE LENGTH FADEN-LÄNGE | CODE ARTIKEL |
|--------|-----|------------------------------|-----------------|
| 1.5 | 4/0 | 45cm | PR532 |
| 2.0 | 3/0 | 45cm | PR533 |

www.medipac.gr

MEDIPAC SA
Industrial Area Kilkis
P.O. BOX 1, GR 61100
Greece

Tel.:+30 2341 071991
Fax:+30 2341 071979
info@medipac.gr

surgical
sutures

PTFE SUTURES - PROFIMED

The pure PTFE monofilament suture

medipac[®]
medical supplies

INVESTING IN
INNOVATION
& HIGH QUALITY