

# Nordent®

Instruments of excellence for the dental industry

BECAUSE WE MAKE THE  
MOST IMPORTANT INSTRUMENTS  
IN THE WORLD...YOURS!

New!  
**XDURA™**

## The Next Generation of Scalers & Curettes

### Cutting Edge Scaler Technology

- New high tech metallurgy advancements
- Uniform particle distribution

### Longest Lasting Scaler in the World

- Spend more time treating patients
- Save money with long lasting scalers

### KEEP EVERYBODY ON THE CUTTING EDGE!

XDURA™ is the hardest, sharpest, longest lasting scaler in the world...period! That means you can spend more time practicing and less time sharpening, you and your patients will be more comfortable throughout the procedure and your office saves money by reducing replacement costs with the longest lasting scaler available. You deserve the best and so do your office and patients...

XDURA™, the next generation of scalers and curettes!

New technology  
XDURA™ stainless  
steel stays sharp  
longer than ANY  
other brand

DuraLite® ROUND  
handle provides  
optimal comfort  
and maximum  
tactile sensitivity

Advanced micro-finishing  
process produces the  
sharpest edge possible

XDURA™ stays sharp  
longer and is easy  
to sharpen!

**XDURA™ SCALERS & CURETTES**

# New! XDURA™

The Next Generation  
of Scalers & Curettes

## A QUANTUM LEAP IN SCALER TECHNOLOGY

- Lasts longer than ANY other brand to reduce your instrument costs
- Consistent, dependable sharpness for better efficiency and less fatigue
- Provide better patient care with sharpness you can trust to perform

## THE NORDENT DIFFERENCE

For more than 40 years, Nordent has been providing innovative products and solutions, designed for real practitioners working in the real world. Our craftsmen take great pride in producing the finest instruments for you, right here in the USA. With Nordent you are assured practical solutions, superior products and a commitment to those who matter most—our customers.

To locate your Nordent representative, review our complete product line, or find practical information for your practice, visit us at [www.nordent.com](http://www.nordent.com)

# Nordent®

NORDENT MANUFACTURING, INC.  
610 BONNIE LANE  
ELK GROVE VILLAGE, IL 60007  
PHONE: 847.437.4780  
FAX: 847.437.4786  
TOLL FREE: 800.966.7336 (US & CANADA)

[www.nordent.com](http://www.nordent.com)




## IT'S WHAT'S ON THE INSIDE THAT COUNTS!

XDURA™ takes advantage of advancements in metallurgy and cryogenic hardening to produce a microstructure not previously possible. All other scalers are made from stainless steels that have large, isolated carbide molecules, making hardness and wear resistance inconsistent. XDURA™ stainless steel features extremely fine carbide particles that are evenly distributed throughout the microstructure to provide you with a blade that is consistently harder, sharper and longer lasting than NY other brand.



## DEPENDABLE SHARPNESS TO THE CORE!



XDURA Scalers & Curettes

Under 500x magnification, these SEM photos clearly show the superior micro-structure and uniform carbide particle distribution of our new XDURA™ Scalers & Curettes.



Leading sharpen-less scaler



All other major brands (440A)

The advanced micro-structure of XDURA™ provides the longest lasting edge possible!



XDURA™ SCALERS & CURETTES